





 **SINMED**



FIXATIONS
CATALOGUE 2021



STABILIZATION OF ENDOTRACHEAL TUBE

IGO-FIX is a disposable, soft band for fixing the endotracheal tube. It provides safe stabilization by preventing accidental extubation. It is made of soft foam covered with a delicate, anti-decubitus material for the skin. The combination of safety and speed of use of the product creates optimal respiratory protection for the intubated patient. It can be used successfully with patients intubated both in the long-term and in an emergency. These bands do not require the use of adhesives, they can be used successfully, with patients with a beard or damaged skin.

IGO-FIX One-piece endotracheal tube holder

Product features:

- Fixation of endotracheal tubes, laryngeal tubes and laryngeal masks
- One-piece
- Quick and comfortable placement
- Can be cut individually – length adjustment to the patient
- Possibility of changing the position
- Comfortable hook strip fastening
- Hole in the band providing immobilisation of the tube
- Additional self-adhesive hook strip on the endotracheal tube enabling its attachment to the band
- Disposable
- Individually packed



REF	LENGTH
010111	72 cm
010112	85 cm

IGO-FIX Two-piece endotracheal tube holder

Product features:

- Fixation of endotracheal tubes, laryngeal tubes and laryngeal masks
- Two-part
- It consists of two elements: neck pad and a tube stabilizing part
- Quick and comfortable placement
- Can be cut individually – length adjustment to the patient
- Possibility of changing the tube position
- Comfortable hook strip fastening
- Hole in the band providing immobilisation of the tube
- Additional hook strip adhesive and silicone strap for the tube enabling its attachment to the band
- Universal size
- Disposable
- Individually packed



REF	PRODUCT DESCRIPTION	REF	PRODUCT DESCRIPTION	REF	PRODUCT DESCRIPTION
01021	1 neck pad 1 tube stabilizing part	01031	1 neck pad 2 tube stabilizing parts	01041	1 neck pad 3 tube stabilizing parts

ETT HOLDER is a reliable fixation of the endotracheal tube, used mainly in emergency medical services. It significantly facilitates and improves the fixing of the tube after its insertion into the trachea.

ETT HOLDER Endotracheal tube holder

Product features:

- Fast and stable fastening of endotracheal tube
- Universal Velcro fasteners enabling effective fixing of the stabilizer around the head
- A special clamping screw that allows you to mount any size of the endotracheal tube
- Teether preventing traumatization of the patient and damage to the tube
- Additional hole for suction without removing the holder
- Special foam on the inside for additional patient comfort



REF	SIZE
ETH-01	UNIVERSAL



TRACHEOSTOMY TUBE HOLDERS

EVO-FIX is a disposable, soft band for fixing the tracheostomy tube. It provides safe stabilization preventing its displacement. Made of soft foam covered with delicate, anti-decubitus material. Adjustable length of straps allows adjustment to any circumference of the patient's neck. Fastening with Velcro helps to stabilize the tube. The combination of safety and speed of use of the product creates optimal protection for patients with a tracheostomy tube.

EVO-FIX Two-piece tracheostomy tube holder

Product features:

- Two-part
- Made of soft and delicate material
- Possibility of changing the tube position
- Comfortable hook strip fastening
- Adjustable length of the band
- Side fastening of both parts hook strip
- Disposable
- Individually packed



REF	SIZE
02009	ADULT

EVO-FIX One-piece tracheostomy tube holder

Product features:

- One-piece
- Quick and comfortable placement
- Possibility of changing the tube position
- Comfortable hook strip fastening
- Adjustable length of the band
- Available version with an additional ventilation hole in the band that allows the skin to breathe
- Disposable
- Individually packed



REF	SIZE	HOLE
02001	ADULT	NO
02002	CHILDREN	NO
02003	INFANT	NO

REF	SIZE	HOLE
020011	ADULT	YES
020021	CHILDREN	YES
020031	INFANT	YES



TORA-FIX a one-piece band for IBP measuring transducers. The wide cuff allows the transducers to be mounted on the patient's arm. Hook strips fasteners enable changes to the position and position of the band repeatedly. It allows the making measurements, even during the transport of the patient.

TORA-FIX One-piece holder of the IBP transducers

Product features:

- One-piece
- Quick and comfortable placement
- Soft and delicate material
- Adjustable length
- Possibility to protect up to 3 transducers
- Comfortable hook strip fastening
- Disposable
- Individually packed



REF	SIZE
06001	10 cm x 30cm
06000	12 cm x 45 cm

NANO-FIX is a one-piece padded fixation for pulse oximeter sensors. It can be applied to the foot or finger. Hook strip fasteners enable the change of the band position repeatedly.

NANO-FIX Fixation for pulse oximeter sensors

Product features:

- One-piece
- Quick and comfortable placement
- Soft and delicate material
- Adjustable length
- Can be used on a foot or finger
- Comfortable hook strip fastening
- Disposable
- Individually packed



REF	SIZE
08001	4 cm x 17 cm
08002	3 cm x 12 cm
08003	1,8 cm x 10 cm




IV-SPLINT is a one-piece immobilization splint that facilitates insertion and infusion into the peripheral arteries and vessels. Sizes available for different age groups to work in every situation. The possibility to adjust the shape to the place of application. Soft and delicate material ensures comfort and safety of the patient.

IV-SPLINT Immobilisation splint

Product features:

- One-piece
- Quick and comfortable placement
- Soft and delicate material
- Inside is a aluminum plate enabling the desired shape
- Elastic straps for stabilization at the site of insertion
- Comfortable hook strip fastening
- Disposable
- Individually packed



REF	SIZE
04001	3 cm x 9 cm
04002	4 cm x 13 cm
04003	5 cm x 16 cm


REF	SIZE
04004	9 cm x 22 cm
04005	9 cm x 33 cm

LIMB-HOLDER is a one-piece band for immobilizing limbs. It prevents the patient from manipulating the punctures and connected wires. Sizes available for different age groups to work in every situation. The integrated handle and a long strap allow fasten of the strap ends to the bedside. Soft and delicate material ensures comfort and safety of the patient.

LIMB-HOLDER Immobilizer

Product features:

- One-piece
- Quick and comfortable placement
- Soft and delicate material
- Adjustable length
- An additional 150 cm long strap with an integrated grip for its attaching to the bedside
- Can be placed on the ankle or the wrist of the patient
- Comfortable hook strip fastening
- Disposable
- The set contains 2 holders



REF	SIZE	REF	SIZE	REF	SIZE
03011	3 cm x 12 cm	03013	7 cm x 22 cm	03015	7 cm x 30 cm
03012	4 cm x 17 cm	03014	7 cm x 25 cm		

PRO-GLOV is a product designed to facilitate the control of disoriented or altered patients requiring a light restraint. Allow to secure the patient's hands to the bed and prevents self-injuries and equipment manipulation. It offers total protection of hands and wrists. Suitable for people who tend to self-injure themselves

PRO-GLOV Protection glove

Product features:

- Quick and comfortable placement
- Soft and delicate material
- Breathable material
- An additional strap with an integrated grip for its attaching
- Disposable



REF	SIZE
2300	UNIVERSAL



VER-FIX PHOTOTHERAPY EYE PROTECTORS

VER-FIX infant Phototherapy Eye Protectors effectively shields a baby's eyes during phototherapy or anytime light protection is required. They ensure maximum protection of the patient's eyes from potentially harmful rays during phototherapy. Delicate material provides comfort for the patient.

VER-FIX One-piece phototherapy eye protectors

Product features:

- One-piece
- Quick and comfortable placement
- Possibility of changing the position
- Comfortable hook strip fastening
- Ability to adjust the size
- Disposable
- Individually packed

REF	SIZE
07001	S
07002	M
07003	L



VER-FIX Two-piece phototherapy eye protectors

Product features:

- Two-piece
- Quick and comfortable placement
- Possibility of changing the position
- Comfortable hook strip fastening
- Ability to adjust the size
- Disposable
- Individually packed

REF	SIZE
07011	S
07012	M
07013	L



SINMED



ul. Graniczna 32B
44-178 Przyszowice
Poland



tel.: +48 (32) 729 82 36
fax: +48 (32) 739 00 45
tel. kom.: 512 997 426



export@sinmed.pl



www.sinmed.pl